

Discovery DMA 850



Site Preparation Guide

Table of Contents


Table of Contents	2
Ideal Setup	3
System Components.....	4
Instrument Measurements	5
Utility Requirements.....	6
Computer Requirements	7-8
Hardware Requirements.....	7
Software Requirements	8
Accessories	9-16
Gas Cooling Accessory (GCA)	9
Air Chiller System (ACS)	10-12
Nitrogen Purge Cooler (NPC)	13
DMA-RH Accessory	14
Air Filter Regulator Assembly	15
Air Compressor Accessory (ACA)	16
Site Preparation Checklist	17
TA Instrument Offices.....	18

Ideal Setup



IDEAL PLACEMENT AND BENCH MEASUREMENTS

Select a location with a rigid laboratory bench that is level and is in a vibration-free environment.



Floor space for ACS: 4 ft²

Floor space for GCA: 4 ft²

Floor space for DMA-RH: 1.5 ft²

Bench width: 61–183 cm
(24–72 in)

Marble table width: 89 cm (35 in)


Marble table depth: 61 cm (24 in)

Marble table distance from the
wall: 31 cm (12 in) min.

System Components



MAIN SYSTEM COMPONENTS



- A. Instrument
- B. Controller
- C. Power Control Box

Instrument Measurements



DMA 850



Height with furnace OPEN: 68.5 cm (27 in)

Height with furnace CLOSED: 53.5 cm (21 in)

Width: 38 cm (15 in)

Depth: 47 cm (18.5 in)

Weight: 34 kg (75 lbs)



POWER CONTROL BOX


Height: 25 cm (10 in)

Width: 18 cm (7 in)

Depth: 43 cm (17 in)

Weight: 5.6 kg (12.5 lbs)

Power cable length:
200 cm (78 in)



Utility Requirements



POWER

- 100–240 VAC, 47–63 Hz, 1400 W
- Safety ground per local regulation
- Neutral to ground <1 volt



European



NEMA 5-15P



GAS

Furnace purge gases: air, nitrogen, argon, helium

Description	Requirement
Air bearing pressure:	• 410–450 kPa (60–65 psig) • Flow rate: 4 L/min
Cooling gas pressure:	175–830 kPa (20–120 psig)
Conditions	• Must be dry • Must be free from oil and dirt
Other	1/8" polyethylene tubing and fittings are supplied in the accessory kit

Computer Requirements



HARDWARE REQUIREMENTS

Description	Requirement
Processor	• Intel® Core™ i5 8400 or better • 2.8 GHz with 9 MB L2 cache
Memory	≥ 16 GB RAM DDR4 2666 SDRAM
Hard drive	≥ 80 GB free space • 1.5 GB required for Full version of TRIOS • 675 MB required for Lite version of TRIOS (without Online help)
DVD (optional)	≥ 48x CD-ROM or DVD (optional for installing TRIOS)
Screen resolution	Required: 1280 x 1024 with 24-bit colors Recommended: 1920 x 1080 with 24-bit colors
Graphic memory	128 MB
Screen (LCD) size	Required: 19" or greater Recommended: 24" wide screen
Network card	Ethernet 10Base T/100 Base TX
Additional Ethernet card(s)	Necessary if connecting the instrument directly and access is needed to the Corporate LAN.
Ethernet Cabling	10/100BaseTX Ethernet hub/switch. Must be EIA-568B Category 5+ UTP
Client-Server Protocol	DHCP
Image Capture (Camera Option)	DirectX 9.0 or higher
TCP/IP ports used	• TCP: 20010, 20011 • UDP: 5050, 5056

Computer Requirements




SOFTWARE REQUIREMENTS

Item	TRIOS
Operating System	• Windows 10 or 11 Enterprise, Ultimate, & Professional • Home version not supported • ≥ 64-bit version
Internet	Internet connection is strongly recommended for ongoing support after installation
Service Pack	Microsoft Operating System Service Pack
Updates	Windows Operating System and associated Microsoft updates must be up to date. Windows 10 must be 1709 or later.
Network	<i>A second network card for corporate connection is recommended.</i> <i>TA Instruments is not responsible for resolving issues associated with connections to your corporate network.</i>
Conflicts	<i>TA Instruments is not responsible for resolving hardware/software conflicts created by the addition of third-party hardware or software to the computer.</i>

Gas Cooling Accessory (GCA) V3



GCA V3 MEASUREMENTS



Dewar diameter: 46 cm (18 in)

Height: 112 cm (44 in)

Width: 46 cm (18 in)

Depth: 79 cm (31 in)

Weight EMPTY: 47 kg (104 lbs)

Weight FULL: 88 kg (194 lbs)

Transfer line: 1.8 m (6 ft)



GCA V3 REQUIREMENTS

Requirements



- 100–240 VAC at 0.9 kVa, 47–63 Hz
- Approved for operation on a 20 A branch circuit with protective conductor (ground)



- Bulk low-pressure liquid nitrogen tank to refill the 50 L GCA dewar
- Pressure: 90 kPa gauge (13 psig) relief valve on tank
345 kPa gauge (50 psig) on fill line
- Use low pressure supply tank only. Recommended filling pressure is 140 to 170 kPa gauge (20 to 25 psig)



Circulator



Power



Cooling



Gas



LN₂



Fluid



Light



Hardware



Software



Temp



Lab



Customer

Air Chiller System (ACS)



ACS-2 MEASUREMENTS WITH CHILLER PANEL



Height: 88.5 cm (35 in)

Width: 67 cm (26.5 in)

Depth: 56 cm (22 in)

Weight: 103 kg (227 lbs)



ACS-3 MEASUREMENTS WITH CHILLER PANEL



Height: 112 cm (44 in)

Width: 67 cm (26.5 in)


Depth: 56 cm (22 in)

Weight: 128 kg (282 lbs)


Air Chiller System (ACS)



CHILLER PANELS



DMA Chiller Panel
P/N 986400.901



Combined DHR/DMA Chiller Panel
P/N 405400.901
P/N 986400.901 with DHR upgrade kit installed

Height: 64 cm (25 in)

Width: 56 cm (22 in)

Depth: 38.1 (15 in)

Weight: 986400.901–6.75 kg (15 lbs)
405400.901–7.25 kg (16 lbs)

Air Chiller System (ACS)



AIR CHILLER SYSTEM REQUIREMENTS

Requirements



- 50 Hz: 207–252 VAC (refer to the serial number plate on the rear of the unit. The ACS is line frequency specific.)
- 60 Hz: 216–252 VAC (refer to the serial number plate on the rear of the unit. The ACS is line frequency specific.)
- 8A, 9A, 9.1A, or 11A (refer to the serial number plate on the rear of the unit)
- US sites require an L6-20 single-phase outlet



NEMA L6-20 plug



- Gas:
 - Air or nitrogen
 - Pressure: 6.9 bar (100 psig)
 - Flow rate: 200 SLPM
 - Temperature: 20–30°C
 - Dew point: -40°C (-40°F) **NOTE:** Dew point is specified at operating pressure. Supplying dryer air at a lower dew point will extend continuous operation.



- Lab Environment (must be below 25°C):
 - 21°C–24°C = Ideal
- Leave 20 cm (8 in) of space in the front and back of the ACS for ventilation.
- 4 ft² required floor space for the ACS.



Circulator



Power



Cooling



Gas



LN₂



Fluid



Light



Hardware



Software



Temp



Lab



Customer

Nitrogen Purge Cooler (NPC)



NITROGEN PURGE COOLER MEASUREMENTS



Height: 38 cm (15 in)

Width: 30 cm (12 in)

Depth: 18 cm (7 in)

Weight EMPTY: 3 kg (6 lbs)

Volume: 2.5 L



NITROGEN PURGE COOLER REQUIREMENTS

Requirements



- Nitrogen only
- Pressure: 172–827 KPa (25–120 psig)
- Flow Rate: 30 L/m



- Relative humidity: 5% to 80% RH from 15°C to 31°C, decreasing to 66% RH at 35°C (non-condensing)
- Operating altitude: 2000 meters minimum



Fill the container with LN₂



Circulator



Power



Cooling



Gas



LN₂



Fluid



Light



Hardware



Software



Temp



Lab



Customer

DMA-RH



DMA-RH ACCESSORY MEASUREMENTS



Height: 38 cm (15 in)

Width: 30 cm (12 in)

Depth: 51 cm (20 in)

Weight EMPTY: 30 kg (67 lbs)



DMA-RH ACCESSORY REQUIREMENTS

Requirements



90–265 VAC, 45–65 Hz



- Dewpoint of -40°C or less
- Pressure: 20 psig
- Flow Rate: 3 L/m
- Nitrogen



15–35°C ambient temperature



Requires distilled water. **NOTE:** Distilled water ONLY. Deionized, HPCL, RO, and tap water are not suitable and may damage the unit and void the warranty if used.



Circulator



Power



Cooling



Gas



LN₂



Fluid



Light



Hardware



Software



Temp



Lab



Customer

Air Filter Regulator



AIR FILTER REGULATOR MEASUREMENTS



P/N 986320.901

Height: 24 cm (9.4 in)

Width: 23 cm (9 in)

Weight: 2.7 kg (6 lbs)

Port size: 1/8"

Distance between mounting holes:
7.5 cm (2.95 in)



AIR FILTER REGULATOR REQUIREMENTS

Requirements



- Inlet pressure: 551–1034 KPa (80–150 psig)
- Outlet pressure: 413 KPa (60 psig)
- Must be dry, and free from oil or dirt



Circulator



Power



Cooling



Gas



LN₂



Fluid



Light



Hardware



Software



Temp



Lab



Customer

Air Compressor Accessory



AIR COMPRESSOR MEASUREMENTS



Height: 23 cm (9 in)

Width: 38 cm (15 in)

Depth: 38 cm (15 in)

Weight: 10 kg (22 lbs)



AIR COMPRESSOR REQUIREMENTS

Requirements



Outlet pressure: 413 KPa (60 psig)



120 VAC, 50/60 Hz, 0.3 Kw or TA Instruments Step-Down Transformer (P/N 203570.001) for 220 VAC environments.



Instrument should not be in an enclosed area; heat build-up occurs if air does not flow freely.



Circulator



Power



Cooling



Gas



LN₂



Fluid



Light



Hardware



Software



Temp



Lab



Customer

Site Preparation Checklist



Discovery DMA 850

	Sufficient bench and floor space for instrument, computer, and cooling accessory (if needed)	
	☐ Bench length: 61–183 cm (24–72 in) ☐ Bench depth: 76 cm (30 in) ☐ Marble table length: 89 cm (35 in) ☐ Marble table depth: 61 cm (24 in)	☐ Floor space for cooling accessory: ☐ GCA: 4 ft ² ☐ ACS: 4 ft ² ☐ DMA-RH: 1.5 ft ²
	☐ Neutral to ground < 1 volt ☐ Instrument power: 100–240 VAC, 47–63 Hz, 1440 W ☐ Pneumatics: 80–150 psig	☐ ACS power: ☐ 207–252 VAC, 50 Hz ☐ 216–252 VAC, 60 Hz
	☐ Computer meets all hardware requirements ☐ I have two network cards installed on the computer to be able to connect to both the corporate LAN and the instrument. ☐ Computer meets all software requirements ☐ Customer's IT personnel has provided Administrative privileges on the controller computer ☐ The Customer's IT personnel will be on site the day of installation	
	Accessory (check all that apply): ☐ Air Chiller System with Chiller Panel ☐ Air Filter Regulator Assembly ☐ Air Compressor Accessory ☐ Nitrogen Purge Cooler ☐ Gas Cooling Accessory ☐ DMA-RH Accessory	
	☐ The Customer assumes responsibility for any damage that occurs when the instrument is moved by someone other than a trained TA Instruments Service Representative.	

I hereby acknowledge that all utility requirements have been met per the checklist above and that they will be ready at the agreed time of installation.

If all utility requirements are not met at the agreed time of installation, additional charges may be incurred for a return Service trip.

Customer _____ DD / MM / YYYY

Company _____ City _____ State _____ Country _____

Please send a signed copy of the completed checklist to your local Service representative.

TA Instruments Offices

For information on our latest products, contact information, and more, see our website at:
<http://www.tainstruments.com>.

To find your local TA Instruments office and contact information, visit
<http://www.tainstruments.com/contact/ta-directory/>

TA Instruments – Waters LLC
Corporate Headquarters
159 Lukens Drive
New Castle, DE 19720
USA

Telephone: 302-427-4000
Fax: 302-427-4001
Email: info@tainstruments.com